Burette Clamp

Bürettenhalter

Georg-August-Universität Göttingen





Burette Clamp

The burette clamp is an aid to secure a burette. It makes sure that the burette is always clamped vertically.

The burette clamp is attached to a lab stand.



The burette clamp is opened using index finger and thumb.



The burette is secured by positioning it between the rubber rollers. Subsequently, the clamp is carefully closed and the burette is ready for use. **Left:** A complete apparatus with burette and burette clamp.

Right: Should no burette clamp be available in the lab, the burette can carefully be secured using a support clamp as shown here.

