

Stoppers

Stopfen

Georg-August-Universität Göttingen



Stoppers

Stoppers are used to cap laboratory equipment.

As needed, glass, plastic or rubber stoppers can be used.





Ground-glass stoppers are used to cap glass apparatuses with ground-glass joints and should be secured with a joint clip.



Plastic stoppers are used to cap glass equipment. They can also be used to seal extraction vessels.



Rubber stoppers can be pierced and used e.g. to secure glass tubes.

Furthermore, they can be used to cap test-tubes.